


**ZOLL M2™**

**ZOLL®**



**Empowering You to  
Save More Lives**


# More Trusted Technology, More Value

ZOLL defibrillators are designed to meet all your resuscitation needs. The ZOLL M2™ is no exception. From patient monitoring to resuscitation with high-quality CPR, the ZOLL M2 embodies proven, reliable ZOLL core technology to ensure optimal patient care.



## More Trusted Technology

The ZOLL M2 aligns with ZOLL's portfolio of resuscitation products that promote consistent, high-quality CPR and high-current defibrillation for adults and children.



## RECTILINEAR BIPHASIC DEFIBRILLATION

ZOLL's unique high-current waveform has been shown to be superior to monophasic waveforms for termination of ventricular and atrial fibrillation. It has been proven most effective in treating high-impedance patients when compared to truncated exponential waveforms at the same or higher energy levels.



## PATENTED 40 MSEC PACING WAVEFORM

Unique, constant-current, 40 msec pacing has the highest capture rate at the lowest average current required, assuring efficacy and patient comfort. Validated with over 4,000 patients in more than 16 studies.\*

\*Studies available upon request.



## CPR FEEDBACK

Delivery of unaided CPR is not good enough. ZOLL has long recognized the importance of CPR feedback. Our CPR help is used worldwide helping clinicians deliver high-quality CPR – the cornerstone of successful resuscitation.



## REUSABLE HANDS-FREE SOFT PADDLES

ZOLL's trusted technology replaces external paddles and allows safer shock delivery with Real CPR Help® included.

## More Value

Designed for clinicians, hospital staff, and emergency responders to offer high-quality patient care, the ZOLL M2 incorporates proven, best-in-class technology, including: defibrillation, pacing, cardioversion, and real-time CPR feedback—all for a lower total cost of ownership, empowering you to help save more lives.

## More Support

The ZOLL M2 professional defibrillator helps you deliver high-quality, resuscitation care to adult and pediatric patients by offering an array of life-saving features:

- CPR feedback to ensure delivery of high-quality CPR and successful resuscitation
- Cardioversion and pacing for patients with tachycardia and bradycardia
- AED mode with proven algorithms for both adults and children

## More Scalable

From a simple shock box to the most complex patient management challenges, the ZOLL M2 offers various options for all patient resuscitation needs, including 12 lead\*, NIBP\*, SpO<sub>2</sub>\*, EtCO<sub>2</sub>\*, and Temperature\*.



\*These parameters will be available in future updates.

## ZOLL MEDICAL CORPORATION

An Asahi Kasei Group Company | 269 Mill Road | Chelmsford, MA 01824 | 978-421-9655 | 800-804-4356 | [www.zoll.com](http://www.zoll.com)

For more information on ZOLL M2, visit [info.zoll.com/m2](http://info.zoll.com/m2)

Copyright © 2018 ZOLL Medical Corporation. All rights reserved. Real CPR Help, ZOLL M2, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners.

Printed in U.S.A.  
MCN IP 1807 0232-05

For subsidiary addresses and fax numbers, as well as other global locations, please go to [www.zoll.com/contacts](http://www.zoll.com/contacts).

# ZOLL®